myBlend

myBLEMISHpen

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ENGLISH

DEFY THE FUTURE OF YOUR SKIN

myBlend is committed to revolutionizing skincare through innovation, sensoriality and beauty-tech. Based on the principles of epigenetics, our mission is to enable every woman to reveal her skin's true potential.

myBLEMISHpen

Thanks to the combined power of blue and red Light, myBLEMISHpen delivers a targeted and non-invasive action to visibly and lastingly reduce the appearance of blemishes. Designed for precise application, it can be used on various areas wherever imperfections appear.

- **Blue LED (465nm) :** helps purify the skin and reduce the appearance of imperfections.
- **Red LED (660nm) :** helps soothe redness and supports skin regeneration.
- **Microvibration (11000 vibrations/min) :** stimulates microcirculation.


Results: visibly reduces imperfections for clearer, smoother and more even-looking skin.

SPECIFICATIONS

LED wavelengths	465nm±5nm (Blue), 660±5nm (Red)
Max. output optical power density	Green: 20.5Mw/CM^2, Red: 11.5Mw/CM^2
LEDs	8 Dual-core LEDs
Operating temperature	0°C to 40°C \| 32W°F to 104°F
Storing and transport temperature	-5°C to 55°C \| 23°F to 131°F
Operating humidity	30% to 95% without condensation
Storage and Transport Humidity	20% to 95% without condensation
Charging supply	5V 2A
Charging port	USB-C
Charging cable	USB-C to USB-A
Charging cable length	1 meter / 3.28 feet

Charging adapter power output	2A
Battery	3.7V, 400 mAh
Battery type	Lithium-ion
Charging time	120 minutes / 2 hours
Recommended use	5 times per day for 8 weeks
Treatment time	8 minutes
Dimensions of the device	162.5 x 24.8 x 24.3 mm
Weight of the device	90g / 3.17oz

ACCESSORIES

myBLEMISHpen | Travel Pouch | Charging cable | User Manual | Quick Start Guide

ADVERSE REACTION

If you have any adverse reaction, such as discomfort, skin pigmentation, skin irritation or other concerns, from using this Device, stop using it immediately and consult your Healthcare Provider.
Skin irritation, bruising, itching, burning sensation on contact, swelling, erythema, dizziness, nausea, and headaches are rare but possible side effects. Immediately stop using the device and contact your Healthcare Provider.

CONTRAINDICATIONS

- DO NOT use the device if you have any active cancer or cancerous lesion.
- DO NOT use the device if you are pregnant, planning on becoming pregnant, or breastfeeding.
- DO NOT use the device on people under the age of 18.

• DO NOT use the device if you have severe acne.
• DO NOT use the device over areas with phlebitis, thrombophlebitis, and varicose veins.
• DO NOT use the device if you have bacterial infections and inflammation of the skin as this may cause symptoms to worsen.
• DO NOT use the device if you suffer from any photosensitive disorder (sensitization to light).
• DO NOT use the device to treat skin that is painful to touch, raised (higher than surrounding skin), crusting, or bleeding. Consult a healthcare professional.
• DO NOT use the device on areas of skin without pigmentation or those affected by vitiligo.
• DO NOT use the device if you are taking any medication or drugs that can cause photosensitivity. If you are not sure about your medication contact your Healthcare Provider.
• DO NOT use the device if you suffer from light-induced headaches
• DO NOT use the device if you have a medical history of seizures or epilepsy that was triggered by light.
Note: We highly suggest taking caution when taking any (skin related) damaging medications. If in doubt, consult your Healthcare Provider.

Note: If you are under the care of a Healthcare Provider, consult your Healthcare Provider before using this device.

WARNINGS

• STOP using the device and consult your Healthcare Provider if you experience an adverse reaction from using the device.
• STOP using the device if you experience any discomfort or difficulty breathing or swallowing occurring from application of the device on the neck.
• STOP using the device if you experience symptoms of hypotension (dizziness or light-headedness) while using the device.
• STOP using the device if it is overheating. Turn off the device immediately.
• STOP your treatment if it feels hot or your skin feels uncomfortably warm.
• DO NOT let children play with the device. Cleaning and charging the device shall not be done by children.
• DO NOT use the device to treat rosacea, moles, warts, open sores, cancerous lesions, or any skin condition.
• DO NOT use the device on spots without pigmentation or areas on the skin that have vitiligo.
• DO NOT shine the light from the LEDs directly into your eyes or someone else's eyes, , as this may cause discomfort.
• DO NOT place the device if you have any open wounds.
• DO NOT use the device outdoors to prevent damage to the

product.
• DO NOT sit, lay down or step on the device. This will break the device.
• DO NOT drop or insert any object into any openings of the device.
• DO NOT use the device if it is not working properly, or if it has been dropped or damaged. For warranty and repairs contact customer support.
• DO NOT use or store the device near hot surfaces.
• DO NOT charge the device when it is close to or might get in contact with water.
• DO NOT submerge the device with water.
• DO NOT store the device in extreme heat or cold.
• DO NOT place the device in direct sunlight.
• DO NOT attempt to modify the device in any way. Modifying or opening the device will void your warranty. the device contains no serviceable parts. This device contains batteries that are non-replaceable
• DO NOT perform maintenance or take apart the product. In case of malfunctioning or damage, contact customer support.
• DO NOT use the device if damaged or when the charging cable is damaged.
• DO NOT use the device with non-manufacturer-approved accessories.
• DO NOT use solvents or abrasive cleaners on the device.
• DO NOT use the device in any other way then specified in this manual.
• DO NOT use the device if you have had a recent surgery or other surgical procedures.
• DO NOT use the device over a bacterial infection as this may result in the worsening of symptoms.
• DO NOT use the device on painful areas. Consult your Healthcare Provider when you experience painful area's during use.
• DO NOT use the device on, near, or by children or those with reduced physical and mental abilities. Close supervision is necessary when this device is used, cleaned or stored near children or those with reduced physical or mental abilities.
• DO NOT use the device with any other USB-C cable than provided by myBlend.
• DO NOT use the device with any other topical other than specified by myBlend.
• DO NOT use the device over areas of skin that lack normal sensation. Only use the device on normal, intact, clean and healthy skin.
• DO NOT use the device on swollen, infected or inflamed area's or skin eruptions.
• DO NOT apply excessive pressure when using the device on the anterior neck.
•DO NOT use the device at a place where other devices

could interfere by emitting electromagnetic energy. This device can be negatively affected by interference of other devices.

Note: The device should only be supplied with 5V DC low voltage. If any part of the main adapter or cable is damaged, discontinue use. Contact myBlend customer support (serviceclient@my-blend.com) to arrange a cable replacement.

PRECAUTIONS

• Do not scratch or damage the light-emitting surface.
• It will void your warranty if you do not operate the device in accordance with the instruction in the manual.
• Read all instructions carefully before using this device.
• Store the device according to the protocol worked out in this manual.
• This device is intended to be used by one person; we recommend to not share it with anyone else for reasons of hygiene.
• The device is designed for cosmetic use only, for individuals with good health. Please contact your Healthcare Provider if you are unsure about your health status.
• The battery of the device is not removable by the user. If the device is to be disposed of, do not put into a fire. Dispose of the device per local regulations.
• Inspect the device and do not use the device if you see any signs of damage such as internal battery leaks. Contact our customer support.
• The device should only be supplied with 5V DC low voltage. If any part of the main adapter or cord is damaged, discontinue use. Contact customer support to arrange a replacement.
• Consult your Healthcare Provider before using the device if you have a fever, acute inflammation, a severe circulatory disorder, light sensitivity, anaesthetic, or desensitized areas.
• We strongly recommend that if you are concerned that you are taking any medication that causes light sensitivity. Consult your healthcare provider.

TRAVEL AND INTERNATIONAL USE

The device's lightweight design makes it a perfect travel companion. Here are some travel tips for you:

• If you are travelling and find that the power outlet does not meet the required specifications, do not charge this device. Device Charging Power Supply Requirements, Input: 100-240vac. at 50/60Hz, Output: 5vdc, 1A
• Confirm with the steward that the Device can be carried and/or used on the airplane prior to flight.

DISCLAIMER

This Device is not intended to diagnose, cure. or prevent any disease. myBlend makes no claims, representations, or warranties regarding the ability of this Device to cure any physical, skin, or mental condition. A qualified Healthcare Provider should be consulted regarding an condition requiring medical attention.

POTENTIAL RISKS

Immediately after treatment skin may appear red. Skin should return to normal within 1-2 hours.
Temporary headaches immediately following the light treatment. Headaches can be treated with an appropriate pain killer and rehydration.

HOW TO USE myBLEMISHpen

GENERAL USE

1. Prior to use, make sure the device is fully charged.
2. Clean the area of skin to be treated.
3. Make sure the device is unplugged and place the device on the skin.
4. Press and hold the power button to turn on the device, LED will turn on.

4.1. In case of low battery, the light next to the charging port will start blinking.
4.2. To control the motor vibration feature, press the vibration button briefly. Each press activates or deactivates will deactivate the vibration, allowing you to adjust the session to your preferences.
4.3. Leave the device on each area for 2 minutes. Every 2 minutes, a stronger vibration will signal that it's time to move to the next area to treat.
4.4. Each session allows treatment of up to 4 areas. The device switches off after a 8-minute treatment automatically. You may stop the treatment at any time by pressing the power button.

4. POWER BUTTON

4.2 VIBRATION BUTTON

5. Take the device off the treatment area.
Clean the Device using a water-based wipe or a damp clean
6. cloth. Ensure the Device is NOT connected to the power supply when cleaning.
7. Repeat these steps based upon the treatment protocol described in the next part.

TREATMENT PROTOCOL

Use the device 5 times per week for a minimum of 4 weeks. After that you can continue using the device as required. Based upon your lifestyle and skin condition. It is not recommended to use the device more than 7 days per week.

CHARGING

• Charge myBLEMISHpen fully for 2 hours before your first use. Only use a 5V adapter or computer USB port (avoid fast chargers as they can damage the battery).
• To charge the device insert the USB-C cable into the device. Insert the USB-A side of the cable into a USB power supply. The LED indicator on the device will start to flash slow or blink slow to indicate charging. When the device is fully charged the LED indicator will turn off.
• Before charging, make sure that the plug and socket are completely dry. Failure to do so may result in electric shock, short circuit, or fire. Note: Always unplug the device before use, the device will not work while being charged.

CARE AND STORAGE

• The lifetime of the Lithium battery can be maximized by allowing it to fully charge and then fully discharge during regular use.
• Store the device at room temperature in dry conditions.
• Store the device in the supplied case or the packaging.
• Make sure the device will not bend, stretch or crush while being stored or when travelled with.
• Avoid storing the device in dusty conditions, in direct sunlight or in damp and humid conditions.
• This device contains no serviceable parts.
• Before each use, with the device switched off, clean the Device following the cleaning instructions below.
IMPORTANT

• The product has a limited lifecycle of two years under nominal use. Following the treatment protocol for two years is considered as nominal use.
• The device must be handled with care.

CLEANING

Clean the Device with a water-based wipe or damp clean cloth. Stubborn stains can be removed by wiping gently with 70% alcohol.
• DO NOT use any other cleansing Devices as these may damage the Device.
• DO NOT clean the device while it is turned on.
• DO NOT clean the device while charging.

TROUBLESHOOTING

If the actions suggested fail to fix the problem, or the problem reoccurs, and your device is within the warranty period please return to myBlend.
• DO NOT attempt to investigate any further.
• DO NOT attempt to modify the device in any way. The device has NO parts you can repair yourself.

Problem	Cause	Solution
The device doesn't switch on.	The device is not fully charged, because the charging cable was not properly inserted into the power outlet.	Insert the charging cable properly into the power outlet.
	The device is not fully charged, because the charging cable was not properly inserted into the device.	Insert the charging cable properly into the device and charge the device.
	The power outlet didn't work during charging.	Try a different power outlet.
	The ON/OFF button has not been pushed one full second.	Push the on/off button for 1 second straight.
The device switches off unexpectedly.	The battery level is too low.	Charge the device.

DISPOSAL

DO NOT dispose of the device with household waste. The Device should be recycled. The Device contains a Lithium-ion battery. For protection of the environment, please safely dispose of the device following your local guidelines.

WARRANTY

The device is warranted to be free from defects in material and workmanship for 24 months from its original purchase date under normal use. This warranty extends only to the original retail purchaser with original proof of purchase and only when purchased from an authorized retailer or reseller. This warranty only extends if the device is used in conjunction with authorized Device accessories. If the Device should become defective within the warranty period, contact the retailer that you purchased the device. This warranty does not cover Devices damaged by:

- Accident, misuse, abuse, or alteration.
- Use with unauthorized accessories or use other than as instructed.
- Connecting to incorrect current and voltage.

MyBlend shall not be responsible for any incidental, special or consequential damages resulting from the use of this Device. All implied warranties, including but not limited to implied warranties of fitness and merchantability are limited in duration to two years from date of original purchase.
This warranty gives you specific legal rights, and you may also have other rights that vary from country to country. Some countries do not allow the exclusion or limitation of incidental or consequential damages, or limitations on how long an implied warranty lasts, so some of the above limitations and exclusions may not apply to you.
This warranty is given in lieu of all other warranties, written or oral, whether expressed by affirmation, promise, description, drawing, model or sample. Any and all warranties other than this one, whether expressed or implied, including implied warranties of merchantability and fitness for a particular purpose are hereby disclaimed.
If the Device is defective within the warranty period please contact the Customer Service Team.
To register your Device and be kept up to date with new and exciting developments from the Customer Service team, please take a moment to sign up on our website www.myblend.com. To complete your registration online you will need to enter the Device batch number. *

SAFETY

The Device has been tested to or complies with the following standards:

IEC/EC 62133:2017+AMD1 2021 Secondary cells and batteries containing alkaline or other non-acid electrolytes - Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications - Part 2: Lithium systems.

EN 55014-1:2006+A2:2011 Electromagnetic compatibility. Requirements for household appliances, electric tools and similar apparatus.

EN 55014-2:2015 Electromagnetic compatibility. Requirements for household appliances, electric tools and similar apparatus. Immunity. Device family standard.

Restriction of Hazardous Substances (RoHS) Directive (2011/65/EU)

Registration, Evaluation, Authorisation and Restriction of chemicals (REACH) directive (EC 1907/2006)

EN 60335-1: 2012+A11: 2014+A13: 2017+A1: 2019 + A14:2019+A2: 2019 + A15: 2021+A16:2023 – Household and similar electrical appliances - Safety - Part 1: General requirements.

EN 60335-2-27:2013+A1:2020+A2:2020– Household and similar electrical appliances - Safety - Part 2-27: Particular requirements for appliances for skin exposure to optical radiation.

EN 62233:2008 – Measurement methods for electromagnetic fields of household appliances and similar apparatus with regard to human exposure.

FAQ

For more information, please visit the FAQ section of our website:

SYMBOLS

Indicates the Device is CE Marked and authorized to sell in Europe. (Lower-Voltage-Directive)

The Waste of Electrical and Electronic Equipment (WEEE), Indicates the Device may not be discarded as household waste (2012/19/EU).

Restriction of hazardous substances.

This warning symbol on the device label is a warning to read the enclosed documents before use and take all safety precautions required.

This symbol on the device label means: Read the user manual before you start using the device.

This symbol means made in China.

This symbol means the upper and lower operating temperatures that the device can be safely used at.

This symbol means the upper and lower operating humidity that the device can be safely used at.

This symbol means that the cardboard packaging that the device came in can be recycled.

CONTACT INFORMATION

Distributed by:
My Blend
serviceclient@my-blend.com
9 Rue du Commandant Pilot,
92200 Neuilly-Sur-Seine

Manufacturer:
Light Tree Ventures Europe B.V.
Saturnusstraat 95, 2516 AG
The Hague, The Netherlands
Web: www.lighttreeventures.com

Date of issuance: February 2026
Version: 1

myblend-HD-85-manual-V1.9.indd 31 3/17/2026 8:56:03 PM

myBlend